



# Fairview Diagnostic Laboratories

University of MN Medical Center University

Fax Supply Order to: 612-273-0183

Lab: 612-273-7838

Clinic Name: Southwest High School Clinic UA4028

Address/Location: 3414 West 47<sup>th</sup> Street, Room E21

Minneapolis, MN

ZIP 55410

Phone: 612-668-3040

Fax: 612-668-3078

Contact Person: \_\_\_\_\_

Date: \_\_\_\_\_

## Supply Request Form

To comply with federal regulation, these supplies are to be used only to collect and transport the specimens sent to Fairview Diagnostic Laboratories for testing.

Qty	Description	Qty	Description
<b>Specimen Collection Tubes/Containers/Supplies</b>			
	Chlamydia/GC PCR <b>swab</b> collection kits (each)		VAC Tube, Purple (EDTA) 3 mL (50 / pack)
	Chlamydia/GC PCR <b>urine</b> collection kits (each)		VAC Tube, Green (PST gel/lithium/hep) 3 mL (50 / pack)
	Fetal Fibronectin Kit (each)		VAC Tube, Lt Blue(3.2% Na citrate) 1.8 mL (100 / pack)
			VAC Tube, Lt Blue(3.2% Na citrate) 2.7 mL (100 / pack)
	Microtainer Tube, Brick Red Plain (TM) 0.6 mL (each)		VAC Tube, Red (no gel) 4 mL clot activator (50 / pack)
	Microtainer Tube, Purple (EDTA) 0.6 mL (each)		VAC Tube, Red (no gel) 6 mL clot activator (50 / pack)
	Microtainer Tube, Green (TMLHG) 0.6 mL (each)		VAC Tube, Red / Gold ring (Gel Serum Separator) 5 mL (50 / pack)
	Microtainer Tube, Red Gel (TMG) 0.6 mL (each)		Vacutainer, 21 GA 1 ¼" Needle (100 / box)
	Stool Hemocult Cards (each)		Venipuncture Needle Pro Holders (HUBS) (100 / box)
	Urine Container (24 Hour w/ collection hat) (each)		Vacutainer Holders, Plain (each)
	Urinalysis Sterile Vacutainer (each)		Vacutainer Brand Luer Adapters (each)
	VAC Tube, Yellow (ACD) 10 mL (each)		
<b>Pathology Supplies</b>		<b>Cytology Supplies</b>	
	Formalin Containers, Prefilled (20ml) (each)		PAP Vials (SurePath) and Cervical Brooms
	Formalin Containers, Prefilled (60 ml) (each)		PAP Cervical Brooms 
	Formalin Containers, Prefilled (90 ml) (each)		PAP Endocervical Brush 
	Specimen Container, Empty, 1 pint size (16 oz.) (ea)		Cytology Spray Fixative
	Specimen Container, Empty, 1 qt size (32 oz.) (ea)		Microscope Slides
	Specimen Container, Empty, 5.25 quart size (168 oz.)		
<b>Culture Supplies (each)</b>			
	Blood/Fungal Culture Bottles, Aerobic, Adult		Stool Culture (Routine) Para-Pak C& S (Orange)
	Blood Culture Bottles, Aerobic, Pediatric		Stool (Giardia & Cryptosporidium) Para-Pak Formalin (Pink)
	Blood Culture Bottles, Anaerobic, Pediatric/Adult		Stool Pak ULTRA ECOFIX for Ova & Parasites (Green)
	Bordetella pertussis Culture (BG Transport Media)		Sterile Plastic Container (each)
	Fungal Isolator tubes (for HIV type patients – Blast Histo Coccid)		Tube, Port-a Cul Anaerobic Transport
	GC Kit- (Thayer Martin plate + CO <sup>2</sup> ampule in bag)		Urine Culture Transport Stabilizer Tube
	NP Wires (Amies with Charcoal Swabs)		Universal Transport Media (comes with swab)
	Swab Culturette (dual swab)		
<b>Miscellaneous Testing Supplies</b>			
	Glucose Tolerance Glucola (50 gram) (# of bottles)		Glucose Tolerance Glucola (100 gram) (# of bottles)
	Glucose Tolerance Glucola (75 gram) (# of bottles)		
<b>Specimen Transport Supplies</b>			
	Bag, Specimen Transport, 6x9 w/ pouch (100 / bundle)		Transport Vial Tubes w/ Snap Cap (AccuSeal)(1000 / box)
	Slide Holders, Plastic, holds 5 slides		Slide Holder, Cardboard, holds 2 slides
<b>Paper Supplies</b>			
	Laboratory Request Form		Pathology Request Form
	Delivery Labels		PAP Exam Request Form
<b>Miscellaneous Supplies* Needed</b>			

\*Additional supplies that are not considered customary laboratory supplies required to collect and transport samples to Fairview Diagnostic Laboratories may have a fee associated with them.

04/24/14

Supply order filled by: \_\_\_\_\_ Date: \_\_\_\_\_

Questions? Please contact Client Services at 612-273-7838